





Filename: 2024-12-04 - 10.25.21.D5000


Sample Info

Well	Conc. [ng/ul]	Sample Description	Alert	Observations
A1	62.0	Ladder		Caution! Expired ScreenTape device (used after two weeks of first use); Ladder
B1		13		Caution! Expired ScreenTape device (used after two weeks of first use)
C1	0.241	14		Caution! Expired ScreenTape device (used after two weeks of first use); Peak out of Sizing Range



Default image (Contrast 100%), Image is Scaled to Sample

A1: Ladder




Sample Table

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
A1	62.0	Ladder		Caution! Expired ScreenTape device (used after two weeks of first use); Ladder

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	7.54	-	773	-		Lower Marker
100	7.01	-	108	11.31		
250	7.93	-	48.8	12.80		
400	8.10	-	31.2	13.07		
600	8.50	-	21.8	13.71		
1000	8.50	-	13.1	13.71		
1500	7.96	-	8.17	12.85		
2500	6.33	-	3.89	10.21		
3500	4.68	-	2.06	7.55		
5000	2.97	-	0.913	4.79		
10000	3.25	3.25	0.500	-		Upper Marker

B1: 13




Sample Table

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
B1		13		Caution! Expired ScreenTape device (used after two weeks of first use)

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
15	7.22	-	741	-		Lower Marker
10000	3.25	3.25	0.500	-		Upper Marker

C1: 14



Sample Table

Well	Conc. [ng/μl]	Sample Description	Alert	Observations
C1	0.241	14		Caution! Expired ScreenTape device (used after two weeks of first use); Peak out of Sizing Range

Peak Table

Size [bp]	Calibrated Conc. [ng/μl]	Assigned Conc. [ng/μl]	Peak Molarity [nmol/l]	% Integrated Area	Peak Comment	Observations
7	0.0290	-	6.35	12.02		Peak outside of Sizing Range
15	6.54	-	671	-		Lower Marker
311	0.212	-	1.05	87.98		
10000	3.25	3.25	0.500	-		Upper Marker